The Pure-Vu[®] Colon Cleansing System Reduces Lifetime Costs and Incidence of Colorectal Cancer (CRC) – A Cost Minimization Analysis

Authors: H. Ephrath¹, J. Voigt, M. Gralnek I.²

1. Motus GI Medical Technologies Ltd, Israel, 2. Emek Medical Center, Afula, Israel

BACKGROUND:

Bowel prep is a major barrier to patient compliance with colonoscopy. Up to 38% of patients do not complete their bowel prep due to poor palatability and/or intolerance of prep volume. Therefore, up to 25% of patients arrive for their colonoscopy with an inadequate prep. This prolongs procedure time, decreases cecal intubation rate, reduces adenoma and neoplasm detection, increases adverse events and requires repeated colonoscopy exams all leading to increased costs. The Pure-Vu[®] System connects to standard commercially available colonoscopes to facilitate intra-procedural cleaning of poorly prepped colons by simultaneously irrigating and suctioning to evacuate feces and fluids.^[1,2,3]

The Pure-Vu System: Fits over most standard colonoscopes, allowing physicians to cleanse the colon in a safe and effective manner to gain clear visualization of the colon mucosa.



AIM:

We performed a cost-minimization analysis simulating the average lifetime costs and incidence of new CRCs comparing colonoscopy using the Pure-Vu System (PVC) vs. standard colonoscopy (SC).

METHODS:

A Markov model was developed comparing the use of PVC in up to 25% of all colonoscopy procedures vs. SC. The model analyzed patients at both average and high-risk for CRC. Patients were cycled through the model based on the probability of finding an adenoma. Probability of follow up colonoscopy was based on colonoscopy findings and the probability of developing CRC based on follow up care. We assumed an inadequate bowel prep rate of up to 25% for the SC cohort and 5% when using PVC. Cost inputs were based on 2017 Medicare reimbursement. Cost of Pure-Vu was set at \$350. Sensitivity analyses over a wide range were performed to identify key drivers of costs and new CRCs.






Table 1: Input Variables Used in the Cost-Minimization Model

Variable	Value (Sensitivity range evaluated)
Compliance screening CRC ⁴	Mean 60%; range (0-100%)
Probability of inadequate prep SOC ⁵	25% (range 5-30%)
Probability of complications colonoscopy ⁶	0.4% (range: 0-1%)
Compliance with surveillance post initial negative screen ⁷	Mean 58%; range (40-90%)
Compliance with surveillance post adenoma identification ⁷	90-95%
Probability of adenoma ⁸	30% (0-50%)
Probability adenoma cancerous ⁹	0.4% (0-1%)
Probability early vs. late stage cancer ¹⁰	85%/15%
CRC – probability over life of patient – average risk ⁴	4.6%
CRC – probability over life of patient – high risk 11	6-6.2%
Life expectancy after diagnosis of early stage CRC ¹²	24 years
Life expectancy after diagnosis of late stage CRC ¹²	3.8 years
Cost colonoscopy no polyp ¹³	\$1,020
Cost colonoscopy polyp ¹³	\$1,315
Cost treatment early stage cancer – year 1^{14}	\$52,640
Cost treatment early stage cancer - ongoing ¹⁴	\$4,425
Cost treatment early stage cancer – last year of life ¹⁴	\$16,840
Cost treatment late stage cancer – year 1^{15}	\$80,640
Cost treatment late stage cancer - ongoing ¹⁵	\$5,960
Cost treatment late stage cancer – last year of life ¹⁵	\$166,670
Cost complications ⁶	\$10,425

Before and After The Pure-Vu[™] System Use








RESULTS:

See Table 2 for overall model outcomes. Using a 60% adherence rate of CRC screening with colonoscopy, the amount of developed cancers in high risk population was reduced from 0.53 to 0.34 comparing the use of SC and PVC, respectively. Under the same adherence rate, PVC was found to be cost saving as compared with SC in both average and high-risk patients (Avg risk patient costs: PVC: \$4,961, SC: \$5,866; High risk patient costs: PVC: \$7,820, SC: \$8,953). Analysis related to repeated colonoscopy procedures due to inadequate prep found that after accounting for a \$350 cost for the device, PVC decreases cost in these procedures by 77-82% with a cost-savings of \$323/\$458 and \$804/\$951 for average and high risk Medicare/Private Payer patients respectively.

Table 2: Assumes 2017 Screening Recommendations of the U.S. Multi-Society Task Force of Colorectal Cancer (MSTF) and 2012 AGA **Recommendations for Surveillance**

	Inadeq. Prep rate	Compliance to colonoscopy	Early CRC – average risk	Late CRC – average risk	Early CRC – high risk	Late CRC – high risk
Including Pure-Vu procedure	5%	60%	2.16%	0.4%	2.9%	0.51%
SOC	25%		3.1%	0.6%	4.6%	0.8%

Graph 1: Average Standard Colonoscopy Cost vs. Average Colonoscopy Cost When Pure-Vu is Utilized (Over Lifetime in High Risk Patients)



CONCLUSIONS:

Improving the inadequate bowel prep rate via Pure-Vu leads to lower costs and a lower projected incidence of CRC in average and high-risk patients over a lifetime time horizon. This is primarily due to earlier identification of adenomas through improved visualization provided by Pure-Vu and polypectomy. Cost savings from Pure-Vu are driven by the reduced need for follow up colonoscopy procedures due to inadequate preps.

- 1. Gross S et al. *Gastrointest Endosc* 2017, doi: 10.1016/j.gie.2017.10.028.
- . Jiménez JP et al. United European Gastroenterol J 2016: 2: Supplement 1. 3. Van Keulen KE et al. United European Gastroenterol J 2017: 5: Supplement 1
- 4. American Cancer Society; Colorectal Cancer Facts and Figures 2017-2019.
- 5. Kingsley J., et al. Cost effectiveness of screening colonoscopy depends on adequate
- bowel preparation rates a modeling study. PLOS One. 2016;11(12):e0167452
- Meester RGS et al. Variation in adenoma detection rate and the lifetime benefits and cost of colorectal cancer screening. A microsimulation model. JAMA 313(23):2349-2358
- 7. Schoen RE, et al. Utilization of surveillance colonoscopy in community practice. Gastroenterology 2010;138(1):73. Doi.1053/j.gastro.2009.09.062
- 3. Kaminski MF, et al. Quality indicators for colonoscopy and the risk of interval cancer. NEJM 2010;362:1795-1803.
- 9. Strum WB. Colorectal adenomas. NEJM. 2016;374:1065-1075.
- 10. Heitman SJ, et al. Prevalence of adenomas and colorectal cancer in average risk individuals: A systematic review and meta-analysis. Clin Gastro Hepatol 2009;7:1272-1278. 15. Mariotto AB et al. Projections of the cost of cancer care in the US: 2010-2020. . Schoen RE, et al. Incidence and mortality of colorectal cancer in individuals with a family history of colorectal cancer. Gastroenterol. 2015;149(6):1438-1445.



12. https://www.johnhancockinsurance.com/life/life-expectancy-tool.aspx 13. Medicare 2017 national average data. Includes: APC5311 + CPT 45378 (diag colonoscopy) + anesthesia

14. Lang K, et al. Lifetime and treatment-phase costs associated with colorectal cancer: Evidence from SEER-Medicare data. Clin Gastro Hepatol. 2009;7:198-204. Use of table 3 data. Costs from year 2006 inflated to 2017 using medical CPI for hospital outpatient. Jrl. Natl Cancer Instit. 2011;103:117-128. Initial costs were inflated to 2017 using the medical CPI for outpatient care.

